

# PREPARING AND RUNNING NORTHERN BLOTTING GEL (HORIZONTAL)

Hongwei Cheng 1/30/02, commented by TCH 1/31/02

## 1. Prepare MOPS/Formaldehyde Agarose Gel

|   |        |
|---|--------|
| Low EEO agarose (Sigma A6013)               | 1.05g  |
| H <sub>2</sub> O                            | 43.4ml |
| Microwave to boil (melt agarose completely) |        |
| 5x MOPS                                     | 14.0ml |
| 37% Formaldehyde                            | 12.6ml |

Cool the gel mix down to approx. 50°C, and pour to an Owls medium size gel box, insert 12-well or 16-well 1.5 cm combs (**Note: you'd better run a two-tiered gel**).

- Mix RNA samples with RNA sample buffer (**see below, ask TCH for your own stock**). (**Note: Usually you load 5-10ug RNA/lane. If you the RNA volume is <10ul, you can use 15ul RNA sample buffer per sample. For RNA vol. >10ul, you may have to use 20ul of RNA buffer**).
- Heat RNA/sample buffer mixture at 65°C for 5 min., then chill tubes on ice until loading.
- [Optional: Load samples and 2ug RNA ladder on the gel].
- Run at 90V for 90 min. with 1x MOPS running buffer.

## Buffer and Solution:

|                 |        |                          |
|-----------------|--------|--------------------------|
| <b>5x MOPS:</b> | 167.5g | MOPS free acid           |
|                 | 27.2g  | NaOAc-3 H <sub>2</sub> O |
|                 | 100ml  | 0.2M EDTA                |
|                 |        | pH to 7.0 with NaOH      |
|                 |        | Q.S. to 4.0 liters       |

**NOTE: It's important to keep the stock solution in dark at 4°C. Discard if it turns yellow.**

## Super RNA sample buffer:

|       |                            |
|-------|----------------------------|
| 2.0ml | 5x MOPS buffer             |
| 3.5ml | 37% Formaldehyde           |
| 10 ml | Formamide                  |
| 4.0ml | 50% glycerol               |
| 31 ul | Ethidium Bromide (10mg/ml) |

---

Total volume = 20 ml. Aliquot and freeze @ -80°C