

Purification of Genomic DNA from Paraffin Sections for PCR (Alkaline Extraction)

Qing Luo 12-7-03; Commented by TCH

Deparaffinization and H & E staining:

- 1) Prepare paraffin-embedded sections (**Note:** we can get nice DNA recovery from one regular slide, i.e. enough DNA for 50-100 PCR reactions).
- 2) Bake the slides at 55-60°C (in an oven or on a heating block) for 1 hour.
- 3) Deparaffinize (or dewax) the slides with xylene (immerse the slides with xylene for 5 min) × 3 times.
- 4) Rehydrate sections by immersing the slides into 100% ethanol for 5 min × 2 >>> 70% ethanol for 5 min × 2 >>> tap water for 5 min. (**Note:** Do not let slides complete dry from this point on).
- 5) Immerse slides in Harris hematoxylin (Fisher) for 2 min (optimization may be needed).
- 6) Rinse slides with water for 1 min × twice (or until water looks clean).
- 7) Immerse slides in 0.3% ammonium hydroxide 10-20 dips.
- 8) Rinse slides with tap water for 1 min × twice.
- 9) Immerse slides in Eosin-Y (Fisher) 10-20 dips (optimization may be needed).
- 10) Wash with water for 1 min × 3 times. (**Note:** You can check the H & E staining under a microscope. Also see the H & E Staining Protocol for other considerations).

Alkaline Extraction of Genomic DNA:

- 11) Remove unwanted tissue from slides with razor blades.
- 12) Add 5µl 0.2 M NaOH/1 mM EDTA to the H & E stained samples, then spread over edge of specimen with razor blade, and use the blades to scrape off the tissue you want, and collect tissue from each slide to a screw-cap tube containing 30µl 0.2 M NaOH/1 mM EDTA. (**Note:** For small areas, you can scrape tissue with 200µl pipet tip instead of razor blade.)
- 13) Vortex the tubes for a few seconds.
- 14) Heat the tubes in a heating block at 85°C for 30 min (or until most of the tissues are dissolved).
- 15) Vortex the tubes for a few seconds.
- 16) Microfuge the tubes at top speed for 2-5 min.
- 17) Transfer 20µl supernatants to 80µl (or 180µl if the sectioned tissue is large) ddH₂O in a new set of tubes. Store DNA at -20°C.

PCR Amplification:

- 18) For PCR amplifications, use 2-3µl the purified genomic DNA as templates for a 25µl reaction. We usually use the touchdown PCR protocol (followed by 30

cycles of regular cycling amplification). (At the end of this stage, you can check 5 μ l of the PCR products on an agarose gel).

- 19) If necessary, the nested PCR amplification can be carried out by using 1 μ l of the first round PCR products (i.e., outside fragments) as templates (for a 25 μ l reaction) and a pair of nest (or inside) primers. Either a touchdown or regular cycling protocol can be used.
- 20) Check 5-10 μ l of the nested PCR products on agarose gels.