

**MANUAL TOUCHDOWN PCR AMPLIFICATION**  
by TCH 8/13/01

PCR Reaction Mix

|                           |               |
|---------------------------|---------------|
|                           | 25 ul Rx      |
| 10x Buffer                | 2.5           |
| dNTPs (10mM each)         | 3.0           |
| DMSO                      | 1.5           |
| Primer #1 (330ng/ul)      | 0.5           |
| Primer #2 (330ng/ul)      | 0.5           |
| ddH <sub>2</sub> O        | q.s. to 25 ul |
| Taq DNA polymerase (BRL). | 0.5           |
| Template DNA (plasmid)    | (10-200 ng)   |
|                           | <hr/> 25 ul   |

PCR Cycling Program (on Hybaid OmnE with "block" control)

|         |  |                |
|---------|--|----------------|
| Stage 1 | 95°C X 2'  | X 1 cycle      |
| Stage 2 | 92°C X 20"<br>64°C X 30"<br>70°C X 30-60" (or 1kb/min) | X 4 cycles     |
| Stage 3 | 92°C X 20"<br>61°C X 30"<br>70°C X 30-60" (or 1kb/min) | X 4 cycles     |
| Stage 4 | 92°C X 20"<br>58°C X 30"<br>70°C X 30-60" (or 1kb/min) | X 4 cycles     |
| Stage 5 | 92°C X 20"<br>55°C X 30"<br>70°C X 30-60" (or 1kb/min) | X 25-35 cycles |
| Stage 6 | 70°C X 5'  | X 1 cycle      |

Load 5-10ul of the PCR product to a 0.8% agarose gel.

**Note:** 1) Lower total cycle numbers are preferred because of lower mutation rate;

2) To obtain larger quantity of DNA, one may set up 2 to 4 reactions (25 ul each);

3) This protocol is suitable for amplification of a specific fragment from a cDNA or genomic library, as well as RT-PCR analysis.